

LAB GLASS DISPOSAL

Date Posted: 9-20-17

BROKEN, & NON-RECYCLABLE GLASS

ACCEPTABLE items include pipettes, test tubes, empty vials, slides, and general lab glass.

Clean and empty bottles and jars.

NO chemicals or chemical residue.

NO plastic caps, needles, blades, light bulbs or hazardous waste.


CONTAINERS - Use a pail or box designated for this purpose. Label it "broken" or "waste" glass". Line the container with a clear trash bag.

DO NOT overfill.

DISPOSAL – Tie the bag and/or seal the box and place next to your trash bin for disposal by custodial staff.



BROWN BOTTLE RECYCLING



Empty 4 L Winchester bottles may be recycled as follows:

STEP 1 – Remove the cap. Caps are not recyclable and must go to the regular trash.

STEP 2 – Bottles must be rid of all chemical residue. Depending on previous contents, this can be done by triple rinsing with an appropriate solvent and/or air-drying overnight in a fume hood.

STEP 3 – Deface the label by marking it with an "X" or the word "EMPTY" so that people who handle the bottle after you know that it has been cleaned.

DISPOSAL - Place next to your recycle or trash bin for pick-up by custodial staff.

